

RELIABLE



Intervertebral fusion device system



Distribuidor exclusivo para Republica Dominicana,
Cuba y toda Centroamérica.

Plaza Las Ramblas, 2do Nivel (Suite 207)
Av. Juan Pablo Duarte, Santiago.
Republica Dominicana
Tel.: +1 (829) 581-4202

e-mail: info@biortesic.com

www.biortesic.com

www.bot.com.do



Distribuidor exclusivo para Republica Dominicana,
Cuba y toda Centroamérica.

Cervical Cage fusion




Product No.	Specifications
11180004	4*14*11
11180005	5*14*11
11180006	6*14*11
11180007	7*14*11
11180008	8*14*11
11180009	9*14*11
11180010	10*14*11
11180011	11*14*11
11180012	12*14*11
11180013	13*14*11


Lumbar Cage fusion




Product No.	Specifications
10180009	32*10*11
10180008	26*10*9
10180007	22*10*7



High resisting dash limit intensity



High resisting limit fatigue intensity
There is no fragments in fatigue experiment



Elastic modulus approach to bony



Screwdriver (Hex)



Trial Implants 5-9



Distractor II



Insertion Device ϕ 3.5



Bone Cartridge



Fixed Screw



Implant Holder



T-Handle (Small)



Mallet



Bone Funnel ϕ 6



9mm Trial Implants



11mm Trial Implants



7mm Trial Implants



6mm Reamer



10mm Reamer



8mm Reamer



Straight Bone Knife



Short Head Curved Bone Rasp



Bone Curette



Short Head Curved Bone Curette (Lower Curved)



Circular Bone Rasp



Bone Rasp (Standard)



Insertion Device ϕ 6



Bone Impactor (Gun Style)